Neuritogenesis and Synaptogenesis Promoted by Acetylcholine-like Motifs in Biodegradable Polymers

Christiane B. Gumera, Jin Gao, Yadong Wang

Department of Biomedical Engineering, School of Chemistry and Biochemistry, and Petit Institute for Bioengineering and Bioscience, Georgia Institute of Technology, Atlanta, GA 30332, USA.

Introduction: Neurotransmitters have been shown to be essential in embryonic and neonatal neuronal development. Acetylcholine (Ach) induces neurite outgrowth and may promote the formation and strengthening of synapses [1]. We hypothesize that a polymer of appropriate Ach content will promote the maximum neurite outgrowth and synapse formation. Here we present a series of new biodegradable polymers that contain varying concentrations of Ach functional groups through the polycondensation of aminoethyl acetate and diglycidyl sebacate (Fig.1). The concentration of the Ach moieties was adjusted by copolymerization with the inert spacer L-leucine ethyl ester.

Fig. 1. Structural similarity between Ach and title polymers.

Methods: Synthesis of monomers: (i) diglycidyl sebacate monomer. The condensation product of allyl alcohol and sebacic acid is oxidized by m-chloroperoxybenzoic acid to produce said monomer. (ii) aminoethylacetate Aminoethylacetate hydrochloride was synthesized according to Liu et al. [2]. Synthesis of Ach-inspired polymers. The following polymers were synthesized in the presence of equimolar $Mg(ClO_4)_2$: $PSA_{100}L_0$, $PSA_{60}L_{40}$, $PSA_{40}L_{60}$, $PSA_{0}L_{100}$ ($A_x = \%Ach$, $L_y = \%$ leucine ethyl ester) Material characterization. The resultant polymers were characterized with FTNMR, FT-IR, DSC, and GPC. Dorsal root ganglia (DRG) culture. Glass coverslips were coated with the corresponding polymer solution (4.62 mg/ml) and dried overnight at 600 mTorr. DRG from 3-day old Sprague-Dawley rats were cultured on the coverslips. Neurite length was analyzed using ImagePro; synaptophysin and neurofilament staining were performed with appropriate antibodies (Biomeda Corp. and Chemicon Intl., respectively).

Results and Discussion: DRG showed significant neurite extension on PSA $_{60}$ L $_{40}$, PSA $_{40}$ L $_{60}$, and PSA $_{0}$ L $_{100}$ (**Fig. 2a**). DRG had few (<5 short neurites) or no neurite growth on PSA $_{100}$ L $_{0}$. More neurites grew on PSA $_{60}$ L $_{40}$ and PSA $_{40}$ L $_{60}$, showed extensive arborization, and were longer than those on PSA $_{0}$ L $_{100}$. Synapse formation was extensive on PSA $_{x}$ L $_{y}$ as indicated by synaptophysin staining (**Fig.2b**). With at least 45 of the longest neurites that were measured on each polymer, the median neurite length on PSA $_{60}$ L $_{40}$ and PSA $_{40}$ L $_{60}$ were comparable and considerably longer than on PSA $_{0}$ L $_{100}$. At day 4, PSA $_{60}$ L $_{40}$ produced neurites up to 2.8 mm (**Fig.3**).




Fig. 2. Ach-mimetic polymers promote neuritogenesis and synaptogenesis. Both DRG above were grown on $PSA_{60}L_{40}$. (a) At day 8, DRG has long, branching neurites extending all around the body. Image obtained by merging multiple 100x photomicrographs. Scale bar = 1mm. (b) Synaptophysin staining of DRG grown on Ach polymer (200x). Numerous synapses formed between neurites, indicated by punctuate red spots. DRG, which express Ach receptors, might recognize the Ach-mimetic functional groups of the polymer.




Fig. 3. Histogram of neurite lengths on DRG cultured for 4 days. DRG on $PSA_{60}L_{40}$ and $PSA_{40}L_{60}$ had median neurite lengths longer than those on PSA_0L_{100} . On $PSA_{60}L_{40}$, however, there were neurites of almost 3 mm compared to those on $PSA_{40}L_{60}$ of 2.4 mm. There were few or no neurites on DRGs grown on $PSA_{100}L_0$. There appeared to be an optimal concentration to stimulate neurite growth. Insufficient Ach led to shorter neurites while over-stimulus with too much resulted in no neurite extension.

Conclusions: We have synthesized a series of polyesters containing various concentrations of Ach-like functional groups. A 60% Ach content appeared to be the most efficient at promoting neuritogenesis and synaptogenesis. We are currently investigating the potential of this new biomaterial in an *in vivo* nerve regeneration model.

References:

- 1. Tata AM. J Neurosci Res 2003;73:227-234.
- 2. Liu J. J Org Chem 2001;66:5655-5663.